Pulse Oximeter

User Manual

(Model: PF-10A)

English

Instructions for Safe Operation


- Make sure that there is no visible damage that may affect user's safety or measurement performance with regard to sensors and clips. It is recommended that the device should be inspected minimally before each use. If there is obvious damage, stop using the device.
- Special attention should be paid while the oximeter is used constantly under the ambient temperature over 37 °C, burning hurt may occur because of over-heating of the sensor at this situation.
- Necessary maintenance must be performed only by qualified service technicians. Users are not permitted to service this device.
- The oximeter must not be used with devices and accessories not specified in User Manual.

Warnings and Cautions

- Explosive hazard—**DO NOT** use the Oximeter in environment with inflammable gas such as some ignitable anesthetic agents.
- **DO NOT** use the Oximeter while the patient is under MRI or CT scanning. This device is NOT MRI Compatible.
- Discomfort or pain may appear if using the Oximeter continuously on the same location for a long time, especially for patient with poor microcirculation. It is recommended that the Oximeter should not be applied to the same location for longer than 2 hours. If any abnormal condition is found, please change the position of Oximeter.
- The light (the infrared light is invisible) emitted from the device is harmful to the eyes. Do not stare at the light.
- The Oximeter is NOT a treatment device.
- Local laws and regulations must be followed when disposing of the device.
- Keep the Oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- Keep this device away from pets, pests or children.
- If the oximeter gets wet, please stop using it and do not resume operation until it is dry and checked for correct operation. When it is carried from a cold environment to a warm and humid environment, please do not use it immediately. Allow at least 15 minutes for Oximeter to reach ambient temperature.
- **DO NOT** operate the button on the front panel with sharp materials or sharp point.
- **DO NOT** use high temperature or high-pressure steam disinfection on the Oximeter. Refer to Chapter 8 for instructions regarding cleaning and disinfection.
- Pay attention to the effects of lint, dust, light (including sunlight), etc.
- Please keep the cable away from children. It can cause strangulation.
- The biocompatibility testing has been performed on the materials in contact with the person in accordance with ISO10993.
- Be placed on the limited training of a lay operator with respect to the ability to intervene and maintain basic safety and essential performance.
- The PATIENT is an intended OPERATOR.
- **Do NOT** dispose of the appliance with the normal household waste at the end of its life, but hand it in at an official collection point for recycling. By doing this you will help to protect the environment.
- Warning about suffocation caused by charging cables.

1 Overview

- 1.1 Intended Use
This Oximeter is intended for measuring the pulse rate and functional oxygen saturation (SpO₂) through a patient's finger. It is applicable for spot-checking SpO₂ and pulse rate of adult and pediatric patients in homes and medical clinics.
- 1.2 Contraindications
No contraindications.
- 1.3 Views

1.4 Features

- Display SpO₂, PR, PI, and Plethysmogram
- Auto power On/Off
- Change between PR and PI
- Over-limit indication and sound
- Mute sound
- Four direction display
- Setting menu (including over-limit setting)
- Pulse beep
- Continuous or spot check measuring mode
- Record list

2 Charging

- Charge the battery before using.
- Connect the device to computer USB or USB charging adapter with USB cable.
- Note:** The device cannot be used during charging, and if choosing a third party charging adaptor (Class II), select one that complies with IEC60601-1 or IEC60950-1.
- Fully charged.
- The filled part represents the remaining power. If the filled part moves from left to right, the device is charging.
- Low battery. Please charge the device
- Note:** Please use the accessories that are original or approved by our company.

3 POWER ON/OFF


- POWER ON:**
Wear the device, it will turn on automatically.
- POWER OFF:**
The device will turn off automatically after 2 seconds.
- On the menu interface, if there is no key operation for about 30 seconds, the device will automatically exit the menu and then shut down.
- On the recording and playback screen, if there is no key operation for 6 seconds, the device will automatically shut down.

4 Start/Stop Measuring


1. Open the clip and put finger inside the clip (make sure the finger is in full contact with the deep inner side of the clip), and then release the clip.
 2. Wait for 2 seconds; the oximeter will power on and start to measure.
 3. The display screen shows the measurement.
 4. Get the finger out, and the device will automatically power off.
- Attentions for measuring:**
- Do not shake the finger and relax during measurement.
 - Do not put wet finger directly into sensor.
 - Avoid placing the device on the same limb which is wrapped with a cuff for blood pressure measurement or during venous infusion.
 - Do not let anything block the emitting light from device, i.e. do not use finger nail polish/paints.
 - Existence of high intensive light sources, such as fluorescence light, ruby lamb, infrared heater or strong sunshine, etc. may cause inaccuracy of measurement result. Please put an opaque cover on the sensor or change the measuring site if necessary.
 - Vigorous exercise and electro-surgical device interference may affect the measuring accuracy.
 - Nail polish may affect the measuring accuracy, and too long fingernail may cause failure of measurement or inaccurate result.
 - If the first reading appears with poor waveform (irregular or not smooth), then the reading is unlikely true, the more stable value is expected by waiting for a while, or a restart is needed when necessary.
 - If the measurements over the limits, there is a reminder sound. You can press the Display key to mute it.

5 Screen

5.1 Indications and Icons

- Icon [25] indicates the counting-down time if the oximeter works at Spot check mode. The total measuring time for Spot check mode is 30 seconds.
- 5.2 Four Directions of the Screen
The oximeter supports to show the screen in four directions. A short pressing of the Display Key can change display direction by 90°, and change PR/PI at the same time. The four display directions are as shown below.

6 Menu Setup

- During measuring, Long pressing Display key can enter the setup menu screen.

SpO ₂ Lo 89	Mode Continuous	Save exit menu
PR Hi 100	Beep On	Restore default
PR Lo 30	Display Always	Setting menu >>>
Setting menu >>>	Setting menu >>>	<<<Setting menu

- Menu operating procedures:
1. Shortly press Display Key to choose the setting item;
2. Long press Display Key to active the setting item, then shortly press it to modify

- the setting parameter;
 - 3. Long press Display Key to confirm the modification and exit from this setting item.
 - 4. Move the setting item to "Exit", and long pressing Display Key to exit from the setup menu.
- Menu settings:**
- **Over-limit settings:** If the SpO₂ or PR value is over the defined limits, the value will flash.
 - **"Beep":** Pulse beep option. If it is set to on, every pulse beat makes a beep.
 - **"Mode":** Set the measuring mode. "Continuous" and "Spot check" for optional, the default is "Spot check".
 - **Spot check mode:** the measuring time lasts 30 seconds with a counting-down indication. The SpO₂ and PR readings will freeze at the end of 30 seconds, the analysis result for the pulse rhythm will be displayed on the screen as well.
 - **Continuous mode:** measurement will start automatically when finger is inserted into the oximeter, SpO₂ and PR readings will be displayed until the finger is removed from the oximeter.
 - **"Display":** The display screen is always on by default. You can set the display to automatically turn off after 5 minutes, 3 minutes, or 1 minute. Wake the screen by pressing the display key.
 - **"Restore default":** Shortly press Display Key to choose "Restore default" and long press Display Key to reset all parameters to their default values.

7 Record List

On power off status, long pressing the Display key shows the record list screen. On record list screen, a short pressing on the Display key can shift the records display, and if there is no key operation for 6 seconds, then the oximeter will power off automatically again.

S: 98 99 98 97			
P: 68 77 75 82			
M1 M2 M3 M4			
Record List			
● A single group of stable readings will be recorded in the record list each time when the oximeter shuts down regardless of spot-check or continuous mode. However, if the time from displaying valid readings to the end of measurement is less than 5 seconds, then no recording will be done.			
● Up to 12 groups of records can be stored in the record list, the newest record is marked as M1, and the oldest record is marked as M12. The new record will override the previous record.			
● When the device is out of battery, it will power off automatically and all readings will be deleted.			

8 Technical Specifications

Classification	Internally powered equipment
The type of protection against electric shock	Internally powered equipment
The degree of protection against electric shock	Type BF applied parts
Electro-magnetic compatibility	Group I, Class B
Environment	
Temperature	Operating: 5 ~ 40°C Storage: -20 ~ 55°C
Relative humidity (non-condensing)	30% ~ 80% 10% ~ 93%
Atmospheric pressure	700 ~ 1060hPa 700 ~ 1060hPa
Degree of dust & water resistance	IP22
Physical	
Dimension	64mm*38mm*28mm
Weight	About 37 g
Display	OLED
Power and supply	
Input	DC 5V ±10%
Battery	Rechargeable Lithium-polymer
Battery life	24 hours for typical use (max)
Charge time	About 3 hours
SpO₂	
SpO ₂ range	Measuring range: 0% ~ 100%
SpO ₂ Accuracy (Arms)	±2% (70% ~ 100%); ±3% (50% ~ 69%); No definition (0% ~ 49%)
Pulse Rate range	30bpm ~ 250 bpm
Pulse Rate accuracy	±2 bpm or ±2%, whichever is greater
SpO ₂ low limit setting range	85% ~ 99% Default setting: 90%
Pulse Rate low limit setting range	30bpm ~ 60 bpm Default setting: 50bpm
Pulse Rate high limit setting range	100bpm ~ 240 bpm Default setting: 120bpm
Sensor	
Wavelength	Dual-wavelength LED sensor with wavelength Red light: 663nm; Infrared light: 890nm
Maximal average optical output power	≤2mW
Ambient light interference	The difference between the SpO ₂ value measured in the condition of indoor natural light and that of darkroom is less than ±1%.
SpO₂ data averaging time	
SpO ₂ data update period	1s
Data update	<10s
Recorded parameters	SpO ₂ , Pulse Rate

9 Maintenance and Cleaning

- 9.1 Maintenance
The expected service life (not a warranty) of this device is 5 years. In order to ensure its long service life, please pay attention to the maintenance.
 - Please change the batteries when the low-voltage indicator lightens.
 - Please clean the surface of the device before using, with 75% alcohol wipes, then let it air dry or wipe it dry. Do not allow liquid to enter the device.
 - If the Oximeter has not been used for more than 7 days, please charge the Oximeter before use.
 - The Oximeter is calibrated in the factory before sale, so there is no need to calibrate it during its life cycle. Any SpO₂ simulators should not be used to validate the accuracy of the Oximeter, they can only be used as functional testers to verify its precision. The SpO₂ accuracy claimed in this manual is supported by the clinical study conducted by inducing hypoxia on healthy, non-smoking, light-to-dark skinned subjects in an independent research laboratory.
- 9.2 Cleaning and Disinfecting Instruction
 - Surface-clean sensor with a soft cloth damped with a solution such as 75% isopropyl alcohol, if low-level disinfection is required, use a mild bleach solution.
 - Then surface-clean with a cloth damped ONLY with clean water and dry with a clean, soft cloth.
- Caution:**
 - High-pressure sterilization cannot be used on the device.
 - Do not immerse the device in liquid.
 - It is recommended that the device should be kept in a dry environment. Humidity may reduce the life of the device, or even damage it.
 - Do not service and maintain while the device is in use.
- 9.3 Cleaning and Disinfecting Instruction
 - Surface-clean sensor with a soft cloth damped with a solution such as 75% isopropyl alcohol, if low-level disinfection is required, use a mild bleach solution.
 - Then surface-clean with a cloth damped ONLY with clean water and dry with a clean, soft cloth.
- Caution:**
 - Do not sterilize by irradiation steam, or ethylene oxide.
 - Do not use the Oximeter if it is damaged.

10 Troubleshooting

Problem	Solution
The SpO ₂ and Pulse Rate value are unstable	Place the finger correctly inside and try again. Keep calm.
Cannot turn on the device	Charge the device.
No display	Charge the device.
"Poor signal" is displayed on the screen. (Only for spot check mode.)	Place the finger correctly inside and try again.
"-." is displayed on the screen.	Place the finger correctly inside and try again. Keep calm.

11 Symbols

Symbol	Description
	Manufacturer
	Date of manufacture
	Importer
	Serial number
	Indicates a medical device that is not to be disposed of as unsorted municipal waste.
	Follow Instructions for Use.
	Type BF Applied Part
	MRI unsafe. Presents hazards in all MR environments as device contains strongly ferromagnetic materials.
	Protected against spraying water and against access to hazardous parts with a tool, per IEC60529
	No alarm system
	CE marking
	Medical device
	Authorized representative in the European community
	UKCA marking
	Authorized Representative in the United Kingdom
	This product complies with the rules and regulations of the Federal Communication Commission.
	Non-ionizing radiation
	This product complies with verpackG.
	Our products and packaging can be recycled, don't throw them away! Find where to drop them off on the www.quefairemedesdechets.fr site (Only applicable for French market).
	Indicator of poor signal, displayed when there is no available reading.

Appendix EMC

The equipment meets the requirements of IEC 60601-1-2:2014.

Guidance and manufacturer's declaration-electromagnetic emission			
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an environment.			
Emissions test	Compliance	Electromagnetic environment-guidance	
RF emissions CISPR 11	Group 1	The Pulse Oximeter uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR 13	Class B	The Pulse Oximeter suitable for use in all establishments, including domestic establishments and those directly network that supplies buildings used for domestic purposes.	
Voltage fluctuations/flicker emissions IEC1000-3-3	N/A	The Pulse Oximeter suitable for use in all establishments, including domestic establishments and those directly network that supplies buildings used for domestic purposes.	
Table 2			
Guidance and manufacturer's declaration-electromagnetic emission			
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge(ESD) IEC61000-4-2	±8 kV contact ±15kV air	≥8 kV contact ±15kV air	Floors should be wood, concrete or ceramic tile, if floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient/burst IEC61000-4-4	±2kV for power Supply lines ±1 kV for input/output lines	N/A	N/A

Surge IEC 61000-4-5	±1kV line (s) to line(s) ±2kV line(s) to earth	N/A	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC61000-4-11	≤5% U _n (≥95% dip in U _i) for 0.5 cycle <40% U _n (60% dip in U _i) for 5 cycles <20% U _n (30% dip in U _i) for 25 cycles <5% U _i (≥95% dip in U _i) for 5 s	N/A	N/A
Power frequency (50Hz/60Hz) magnetic field IEC61000-4-8	3A/m	3A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: U_i is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity			
The Pulse Oximeter is intended for use in the electromagnetic environment specified below. The customer or the user of the Pulse Oximeter should assure that it is used in such an electromagnetic environment.			
Immunity test	IEC60601 test level	Compliance level	Electromagnetic environment-guidance
Portable and mobile RF communications equipment should be used no closer to any part of The Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance d=1.2 √P (P=800MHz to 800MHz) Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).			
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol.			
Conducted RF IEC61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
Radiated RF IEC61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	a: Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, and electromagnetic site survey should be considered. If the measured field strength in the location in which the Pulse Oximeter is used exceeds the applicable RF compliance level above, the Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Pulse Oximeter.
b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

Rated maximum output power of W(watts)	Separation distance according to frequency of transmitter M(Meters)		
	150kHz to 80MHz d=1.2 √P	80MHz to 800MHz d=1.2 √P	80MHz to 2.5GHz d=2.3 √P
0,01	N/A	0.12	0.23
0,1	N/A	0.38	0.73
1	N/A	1.2	2.3
10	N/A	3.8	7.3
100	N/A	12	23

- For transmitters rated at a maximum output power not listed above, the recommended separation distance in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
- NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
- Copyright**
This manual is written by our company and all rights reserved. Without our company's prior written consent, no part of this manual may be reproduced or copied in any form or method.
- Illustration**
All illustrations provided in this manual are for reference only, and the settings or data in the illustrations may not be exactly the same as the actual display you see on the product.

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